



ASPEN IMAGING
HEALTHCARE

Europa85

Portable X-ray System Unit




DESIGNED FOR MAXIMUM PORTABILITY AND STREAMLINED WORKFLOWS—PERFECT FOR CARE PROVIDERS ON THE MOVE

The Europa85 is Aspen Imaging Healthcare's compact yet powerful portable X-ray solution, engineered for mobility and precision. With an ultra-light weight of just 7.9 lb, it's ideal for mobile healthcare teams, and small practice settings that demand flexibility without compromising image quality.

Unlike conventional portable X-ray devices, the Europa85 allows images to be captured and reviewed directly on the device. Once reviewed, these images can be effortlessly transmitted wirelessly to a PC through its integrated workstation**, enhancing workflow and efficiency.

DEVICE FEATURES



* When using Mirror View, the image transmitted via Miracast is not intended for diagnostic purposes.

**Available under specific conditions.



WHAT'S INCLUDED

- Europa85
- Carrying case
- Handswitch
- Charger
- Stylus pen
- Carrying strap



ALL-IN-ONE PC WORKSTATION

Integrated INTEL PC with
Windows 10+ compatible.
preloaded with DR Software*



EXTREME LIGHT WEIGHT

8 lb for ultimate
portability and convenience



REMOTE SUPPORT



Pre-loaded Remote Assistance
Software (AspenView)**

SPECIFICATIONS

Europa85

POWER OUTPUT	425 W
kV / mA	50 ~ 85 kV / 2~5 mA
EXPOSURE TIME	0.03 ~ 2 sec
COLLIMATOR	Multi-leaf collimator
ANODE TYPE	Stationary anode
ANODE ANGLE	16°
FOCAL SPOT	0.8 mm
HEAT CAPACITY	8,000 J
BATTERY CAPACITY	22.2 VDC / 2,200 mAh
OPERATION	LCD touch screen / PC
WEIGHT	About 3.6 kg / 7.9 lb
DIMENSIONS	266 mm (W) x 167 mm (H) x 156 mm (D)

HANDHELD MODE



* Only preloaded when purchasing our panel.

**AspenView is included only when the panel is purchased.



ASPEN IMAGING
HEALTHCARE

801 Jupiter Rd Suite 200, Plano, TX 75074, USA
sales@aspenimaging.com | (+1) 214-257-0113

WWW.ASPEN-IMAGING.COM

DISTRIBUTED BY

About Aspen Imaging

Aspen Imaging Healthcare is a Texas-based Original Equipment Manufacturer (OEM) providing complete imaging solutions for x-ray distributors and custom designed, private-labeled solutions for OEMs across the Americas.

Founded in 2015, Aspen Imaging has become one of the fastest growing Medical Imaging companies with offices in the United States, Mexico and South Korea.

