



ASPEN IMAGING
HEALTHCARE

Europa90

Portable X-ray System Unit




DELIVERS ENHANCED CONTROL, PRACTICAL IMAGE PERFORMANCE, AND SPEED—IDEAL FOR FOCUSED DIAGNOSTIC ENVIRONMENTS

The Europa90 is a versatile solution crafted to effortlessly integrate into any workspace. Showcasing a compact and featherweight design, its wireless and battery-powered system enables operators to capture X-ray images without the need for the large physical space demanded by traditional radiographic systems.

In contrast to standard portable X-ray devices, the Europa90 allows images to be acquired and reviewed directly on the device. Once reviewed, these images can be seamlessly transmitted wirelessly to a PC via its built-in workstation**, streamlining processes and optimizing efficiency.

DEVICE FEATURES



* When using Mirror View, the image transmitted via Miracast is not intended for diagnostic purposes.

**Available under specific conditions.



WHAT'S INCLUDED

- Europa90
- Carrying case
- Handswitch
- Charger
- Stylus pen
- Carrying strap



ALL-IN-ONE PC WORKSTATION

Integrated INTEL PC with Windows 10+ compatible, preloaded with DR Software*



EXTREME LIGHT WEIGHT

10.8 lb for ultimate portability and convenience



REMOTE SUPPORT


Pre-loaded Remote Assistance Software (AspenView)**

SPECIFICATIONS

Europa90

POWER OUTPUT	450 W
kV / mA	50 ~ 90 kV / 2 ~ 5 mA
EXPOSURE TIME	0.04 ~ 1.3 sec
COLLIMATOR	Multi-leaf collimator
ANODE TYPE	Stationary anode
ANODE ANGLE	16°
FOCAL SPOT	0.5 mm
HEAT CAPACITY	10,000 J
BATTERY CAPACITY	22.2 VDC / 2,200 mAh
OPERATION	LCD touch screen / PC
WEIGHT	About 4.9 kg / 10.8 lb
DIMENSIONS	367 mm (W) x 201 mm (H) x 188 mm (D)

HANDHELD MODE



* Only preloaded when purchasing our panel.

**AspenView is included only when the panel is purchased.



ASPEN IMAGING
HEALTHCARE

801 Jupiter Rd Suite 200, Plano, TX 75074, USA
sales@aspenimaging.com | (+1) 214-257-0113

WWW.ASPEN-IMAGING.COM

DISTRIBUTED BY



About Aspen Imaging

Aspen Imaging Healthcare is a Texas-based Original Equipment Manufacturer (OEM) providing complete imaging solutions for x-ray distributors and custom designed, private-labeled solutions for OEMs across the Americas.

Founded in 2015, Aspen Imaging has become one of the fastest growing Medical Imaging companies with offices in the United States, Mexico and South Korea.

